

Declaration of Sterilization

This is to declare that our gloves specified below have been subjected to a sterilization process in compliance with the relevant national and international standards, and the sterilization effect meets the specified sterile requirements for medical/pharmaceutical industrial use.

KINGSURE CSM Material Glove

Can be sterilized in steam autoclave up to 60 cycles under below conditions:

≤121°C, 30mins

≤130°C, 10mins

Can be sterilized with hydrogen peroxide: 300000ppm / 24h

Can be sterilized with peracetic acid: 15% / 24h

Can be sterilized with Ozone: 200ppm / 72h

Can be sterilized with gamma rays

KINGSURE EPDM Material Glove

Can be sterilized in steam autoclave up to 150 cycles under below conditions:

≤121°C, 30mins

≤130°C, 10mins

Can be sterilized with hydrogen peroxide: 300000ppm / 24h

Can be sterilized with peracetic acid: 15% / 24h

Can be sterilized with Ozone: 200ppm / 72h

Can be sterilized with gamma rays

Kingsure Bio-tech (Shanghai) Co., Ltd



2026-02-05